



The smartest solution  
for temperature monitoring  
of your thermo-controlled chambers



Access all temperature data  
from your smartphone, tablet or computer  
through a secured login access

Alarm notification by SMS, email or phone call  
in case of temperature excursions

Measurement monitoring and metrology





[www.jri-maxant.com](http://www.jri-maxant.com)



SiriusWeb allows to monitor many types of laboratory equipments:

refrigerators, freezers, cold rooms, incubators..., and ambient conditions in production or storage areas.

Different types of alarm notification according to the SiriusWeb service package subscribed

Lite	Gold
	  
Measurements and history available over a rolling 30-day period	Measurements and history available over a rolling 12-month period



The **SiriusBox** measures and records temperature, hygrometry or door openings.

Measurements data are forwarded to the SiriusWeb cloud through your web connection.

In case of thresholds excursion,

**SiriusWeb** sends alarm alerts by e-mail, SMS text message and phone call (according to the level of service you have subscribed).


The SiriusBox integrates an internal buzzer for a local alert (GOLD level service only).


**A comprehensive and automated solution** from measuring to alert notification.

**Easy to implement and use**

**Easier metrology services**

ISO17025 calibration can be performed by simple probe exchange. Your metrological documents can be stored on your SiriusWeb space.





4 <sup>TH</sup> GENERATION DIGITAL PROBES	
Temperature	from -55°C to +80°C
Low Temperature	from -200°C to 0°C
Temperature/ Hygrometry	from -45°C to +85°C/ from 0 to 100%RH

OPERATION

FEATURES